

The Los Alamos Immuno-Lymphocyte Proliferation Test

Babetta L. Marrone, Ph.D.
Group Leader, Cytometry Group
Life Sciences Division
Los Alamos National Laboratory
blm@lanl.gov


LANL LPT Program

- **Development of the Immuno-LPT and Application to Medical Surveillance of LANL Beryllium Workers**
 - » **>300 workers**
 - » **Informed consent(s) in place**
 - » **Testing has begun**
 - » **100 workers tested so far**
 - » **Split samples sent to NJC**
 - » **NJC CBD samples sent to LANL for test validation and research**




Goal and Approach

- **Develop LPT with greater predictive value for CBD**
- **Base LPT on Flow Cytometry**
 - **More sensitive**
 - single cell analysis
 - **More specific**
 - lymphocyte subset analysis
 - **More informative**
 - several measures obtained from the same sample




The Immuno-LPT

- **Rationale: CBD is an immune disease, involving T helper cells**
- **Measure CD3, CD4, CD8, and proliferation in same sample using Multiparameter Flow Cytometry**
- **Express % proliferating cells in terms of cell subset**



The LPT Measures Cell Proliferation in Response to Beryllium



Two Immuno-LPTs Have Been Developed



- **Screening Immuno-LPT**

- BrDU incorporation following Be challenge in CD4 and CD8 subsets
- *Purpose:* Analogous to Thymidine incorporation assay

- **Confirmatory Immuno-LPT**

- DNA content assay following Be challenge in CD4 and CD8 subsets
- *Purpose:* Confirm proliferation in BrDU-positives

Average % Replication after 7d Culture with 10 μ M Be and Human Serum (Difference from Media Control)



“Other Sensitives”

<u>#</u>	<u>Exp.*</u>	<u>BrDU</u>		<u>DNA</u>	
		<u>CD4</u>	<u>CD8</u>	<u>CD4</u>	<u>CD8</u>
1)	F	0.2	0	0.4	2.5
2)	F/C	0.1	0	0.2	2.2
3)	F	1.7	2.2	0.2	4.2
4)	F/C	0	0.9	0	4.4
5)	C	0.2	0.9	0.8	4.3
6)	NE	0	0	0	5.0

Assay results from 10 or 100μM Be and human serum

Exp*:

F= former, C=current, NE= nonexposed

Conclusions

- **2 Immuno-LPTs based on flow cytometry have been developed**
 - screening
 - confirmatory
- **Lack of a confirmed CD4 response may distinguish CBD from sensitive**
 - NEED MORE DATA
- **Plasma (autologous) maximizes CD4 cell response**
- **Serum (non-autologous) maximizes CD8 cell response**